

ASKIR C30

SUCTION UNITS

13



MAIN APPLICATIONS

Dental Practice

Emergency Dept.

Endoscopy

Gastroenterology

General Surgery

Obstetrics

CE 0123

ASKIR C30 is a portable suction unit on castors for the aspiration of body liquids, oral, nasal and tracheal aspiration in adults or children. For professional use and with a powerful aspiration of max 40 l/min, ASKIR C30 comes with double collection jars and on a stand with five castors for easy transport. Main unit equipped with vacuum-meter (bar & kPa) and vacuum regulator.

AVAILABLE MODELS and ACCESSORIES INCLUDED

	BASIC	FS
Liquid Collection Jars with overflow valve system (different options, see below)	o	o
8x14 Silicone Tubes (autoclavable) + Conical Connector ø 10-11-12mm	o	o
Antibacterial & Hydrophobic Filter (single-patient)	o	o
Power Cord with Schuko plug	o	o
Footswitch (may be ordered at a later time)		o

TECHNICAL FEATURES

Motor	Oiless and maintenance-free piston pump
Power Feeding	230V-50Hz
ISO 10079-1 Classification	HIGH VACUUM / HIGH FLOW
Max Vacuum (adjustable)	-0.80 Bar -80 kPa -600 mmHg
Max free air flow rate	40 l/min
Noise Level	61,5 dBA
Duty cycle	Non-stop operation
Weight	Kg 6.20
Size	cm 32 x 99 x 30
Years of Warranty	2
Shipping carton	1
Place of Manufacturing	Italy

Italian
Medical
Touch
MADE IN ITALY

AVAILABLE CONFIGURATIONS

RE 410250	ASKIR C30 BASIC	with 2x2000ml autoclavable jars in polycarbonate
RE 410250/01	ASKIR C30 FS	with 2x2000ml autoclavable jars in polycarbonate and footswitch
RE 410250/10	ASKIR C30 BASIC	with 2x2000ml Flovac® disposable liners and reusable containers
	— Flovac® disposable liners and jars in APEC® autoclavable up to 143°C available upon request	