

# NEW HOSPIVAC 400

SUCTION UNITS

15



## MAIN APPLICATIONS

Dental Practice

Emergency Dept.

General Surgery

Gynecology

Liposuction

Neurosurgery

Obstetrics

Operating Theatre

CE 0123

NEW HOSPIVAC 400 has been designed for professional aspiration of bodily fluids, gases, tissues or bones of the patient during or after surgery. The state-of-the-art oilless pump guarantees high performances with excellent suction capacities and max vacuum reachable within a few seconds, with no need for maintenance. A clear dashboard along with a full range of accessories makes it the ideal device for surgical suction and liposuction.

## AVAILABLE MODELS and ACCESSORIES INCLUDED

	BASIC	FS	FULL
Liquid Collection Jars with overflow valve system (different options, see below)	o	o	o
8x14 Silicone Tubes (autoclavable) + Conical Connector ø 10-11-12mm	o	o	o
Antibacterial & Hydrophobic Filter (single-patient)	o	o	o
Power Cord with Schuko plug	o	o	o
Footswitch (may be ordered at a later time)		o	o
Electronic Change-Over System from jar to jar			o

## TECHNICAL FEATURES

Motor	Oilless and maintenance-free piston pump
Power Feeding	230V-50Hz - other configurations on request
ISO 10079-1 Classification	HIGH VACUUM / HIGH FLOW
Max Vacuum (adjustable)	-0.90 Bar -90 kPa -675 mmHg
Max free air flow rate	90 l/min
Noise Level	46,4 dBA
Duty cycle	Non-stop operation
Weight	Kg 20
Size	cm 46 x 85 x 42
Years of Warranty	2
Shipping carton	1
Place of Manufacturing	Italy

Italian  
Medical  
Touch  
MADE IN ITALY

## AVAILABLE CONFIGURATIONS

RE 410350	NEW HOSPIVAC 400 BASIC 2	with 2x2000ml autoclavable jars in polycarbonate
RE 410350/09	NEW HOSPIVAC 400 FS 2	with 2x2000ml autoclavable jars in polycarbonate
RE 410350/01	NEW HOSPIVAC 400 FULL 2	with 2x2000ml autoclavable jars in polycarbonate
RE 410350/05	NEW HOSPIVAC 400 BASIC 4	with 2x4000ml autoclavable jars in polycarbonate
RE 410350/10	NEW HOSPIVAC 400 FS 4	with 2x4000ml autoclavable jars in polycarbonate
RE 410350/18	NEW HOSPIVAC 400 FULL 4	with 2x4000ml autoclavable jars in polycarbonate
RE 410350/08	NEW HOSPIVAC 400 FLOVAC BASIC 2	with 2x2000ml Flovac® disposable liners and reusable containers
RE 410350/03	NEW HOSPIVAC 400 FLOVAC FS 2	with 2x2000ml Flovac® disposable liners and reusable containers
RE 410350/11	NEW HOSPIVAC 400 FLOVAC FULL 2	with 2x2000ml Flovac® disposable liners and reusable containers
RE 410350/27	NEW HOSPIVAC 400 FLOVAC BASIC 3	with 2x3000ml Flovac® disposable liners and reusable containers
RE 410350/28	NEW HOSPIVAC 400 FLOVAC FS 3	with 2x3000ml Flovac® disposable liners and reusable containers
RE 410350/25	NEW HOSPIVAC 400 FLOVAC FULL 3	with 2x3000ml Flovac® disposable liners and reusable containers
	Jars in APEC® autoclavable up to 143°C available upon request	