



For The Long Run®

Feel The WOODWAY Difference.

Maximum Versatility - The PPS Plus

The WOODWAY PPS Plus treadmill is the ideal device for patients who require product flexibility for therapy. The fully adjustable gas assist handrails are infinitely customizable in height and width to provide a security level unmatched by more traditional stationary parallel handrails.

The generous number of features included in the WOODWAY PPS Plus series offers maximum variability providing solid comfort and accurate support for patient and therapist alike.



Img. PPS 55 Plus

WOODWAY User System (WUS) Display

- LCD touch screen operator panel
- 18 pre-set permanent programs
- Medical test cycles, e.g. Conconi, Cooper, Bruce
- 200 individual custom programs
- Pulse control with all programs
- Remote speed, stop and elevation control
- Speed, distance, elevation, heartrate, time and calorie readout



Img. WOODWAY User System (WUS)

Better for your body:

- Joint-protecting slat-belt technology (Point elastic technology for an optimal impact absorption)
- Smooth, slip-free acceleration, zero start with 0.1 km/h speed increments
- Unparalleled running experience

Better for the environment:

- Significantly reduced energy consumption (low-friction ball bearing transportation system)
- Longer treadmill life
- Factory renewal program

Better for your bottom line:

- Considerably reduced maintenance, spare-parts and repair costs (no replacement of deck or running surface required)
- Significantly reduced operational costs
- Usage possible for patients of all sizes (maximum patient weight up to 360 kg)

Physical Specifications

Belt Type	60 individual slats
Running Surface	9,5mm vulcanized rubber (43-47 shore hardness)
Drive System	114 precision ball bearings with 10 guide bearings (4mm lateral tolerance)
Drive Motor	brushless servo (1,5 kW)
Total Dimensions	PPS 43 184 x 94 x 150 cm (LxWxH) PPS 55 184 x 106 x 150 cm (LxWxH) PPS 70 184 x 121 x 150 cm (LxWxH)
Weight	PPS 43 220 kg PPS 55 240 kg PPS 70 270 kg
Power Supply	220-240V 50/60 Hz (16A fuse)

Performance Specifications

User Weight Cap.	180 kg Running 360 kg Walking
Running Surface (usable)	PPS 43 172(157) x 43 cm (LxW) PPS 55 172(157) x 55 cm (LxW) PPS 70 172(157) x 70 cm (LxW)
Speed Range	0 - 24 km/h (0.1 km/h resolution)
Reverse	0 - 10 km/h at 10% incline max
Elevation Range	0 - 25% (0.1 % resolution)

WOODWAY-Warranty	7 year treadmill frame 5 year function of the running belt* 3 year drive and motor as well as all non-frictional parts* 1 year full warranty incl. labor
------------------	---

* (The above warranty is only valid in conjunction with the performance of a WOODWAY-certified annual maintenance program.)

Authorized WOODWAY Distributor:

Please check our valid pricelist for available options.

Full Compliance with the Legal Requirements for Medical Products in the European Economic Area (CE marking)
 WOODWAY PPS Series treadmills comply with the Medical Devices Directive (directive 93/42/EEC, last amended by 2007/47/EC).